# PH2103: Physics Laboratory II

(PH2203: Physics Laboratory III)

#### Instructors:

Bhavtosh Bansal, Bheemalingam Chittari, Partha Mitra Rumi De

# PH2103: Physics Laboratory II

(PH2203: Physics Laboratory III)

#### **Instructors**:

Bhavtosh Bansal, Bheemalingam Chittari, Partha Mitra Rumi De

#### **Optics**

(Bhavtosh Bansal and Partha Mitra)

#### Modern Physics

(Bheemalingam Chittari and Rumi De)

### Experiments:

- 1. Electron Diffraction,
- 2. Photo Electric effect,
- 3. Velocity of Light,
- 4.Frank Hertz,
- 5. Stefan Boltzmann

Lab Notes format for each experiment:

- Title:
- Aims(s):
- Working Principle/Formula:
- Experimental Setup/Schematic Diagrams:
- Data/Readings/Table/Plot:
- Analysis:
- Source of Errors:
- Conclusions:

# PH2103: Physics Laboratory II

### Modern Physics Lab

#### Source of Errors:

- Systematic Error Instrument resolution,
- Environmental factors Temperature, secondary source feed back, etc..
- Not sufficient readings Fluctuations in readings/measurements, etc..
- Human Error
- Faulty Instruments

1. Electron Diffraction,

#### 1. Electron Diffraction,

Particles have wave properties in addition to their familiar particle properties.

Louis de Broglie in 1924

$$\lambda = \frac{h}{P}$$

 $\lambda$ : (Wavelength), h: (Plank'sConstant) and P: (Momentum)

#### 1. Electron Diffraction,

Particles have wave properties in addition to their familiar particle properties.

Louis de Broglie in 1924


$$\lambda = \frac{h}{P}$$

 $\lambda$ : (Wavelength), h: (Plank'sConstant) and P: (Momentum)

This particle wave nature is confirmed from the experimental observation on the <u>diffraction of electrons</u> in crystalline Nickel structure

#### 1. Electron Diffraction,


The regular arrangement of atoms in a single crystal can be understood as an array of lattice elements on parallel lattice planes.



Lattice plane spacings in graphite

#### 1. Electron Diffraction,


The regular arrangement of atoms in a single crystal can be understood as an array of lattice elements on parallel lattice planes.




Lattice plane spacings in graphite

When we expose such a crystal lattice to *monochromatic x-rays* or *mono-energetic electrons*, and, additionally assuming that those have a wave nature, then each element in a lattice plane acts as a "scattering point", at which a spherical wave- let forms.

#### 1. Electron Diffraction,






- These spherical wavelets create a "reflected" wave front. The wavelength λ remains unchanged with respect to the "incident" wave front, and the radiation directions which are perpendicular to the two wave fronts fulfil the condition "angle of incidence = angle of reflection"
- The constructive interference arises in the neighbouring rays reflected at individual lattice planes when their path differences  $\Delta = \Delta_1 + \Delta_2 = 2d\sin\vartheta$


$$n\lambda = 2d\sin\theta$$

n = 1,2,3... d =lattice plane spacing,  $\vartheta =$ diffraction angle

#### 1. Electron Diffraction,

Diffraction is observed in two ways:






Davisson and Germer in 1927

# $d_1$ $d_2$


#### 1. Electron Diffraction,

Diffraction is observed in two ways:



#### 1. Electron Diffraction,

Diffraction is observed in two ways:



#### 1. Electron Diffraction,





Fig. 4: Schematic sketch for determining the diffraction angle. L=13.5 cm (distance between graphite foil and screen), D: diameter of a diffraction ring observed on the screen  $\vartheta$ : diffraction angle For meaning of F<sub>1</sub>, F<sub>2</sub>, C, X and A see Fig. 5.






Fig. 5: Experimental setup (wiring diagram) for observing the electron diffraction on graphite. Pin connection:

 $F_1$ ,  $F_2$ : sockets for cathode heating


C: cathode cap


X: focusing electrode

A: anode (with polycrystalline graphite foil see Fig. 4)

#### 1. Electron Diffraction,


- The electrons emitted by the hot cathode a small beam is singled out through a pin diagram.
- After passing through a focusing electron-optical system the electrons are incident as sharply limited monochromatic beam on a polycrystalline graphite foil.
- The atoms of the graphite with different space lattice which acts as a diffracting grating for the electrons.
- On the fluorescent screen appears a diffraction pattern of two concentric rings which are centred around the indiffracted electron beam







#### 1. Electron Diffraction,






Since the microcrystals are in all possibe orientations, the diffraction pattern consists of concentric cones with diffraction angles  $2\theta$ . These can be recorded with a cylindrical film or area detector. The intensities are measured as a function of  $\theta$ .








# $d_1$ $d_2$

#### 1. Electron Diffraction,



$$\lambda = \frac{h}{P}$$


• The diameter of the concentric rings changes with the wavelength  $\lambda$  and thus with the accelerating voltage U

e. U = 
$$\frac{1}{2}mv^2 = \frac{1}{2m}(mv)^2 = \frac{1}{2m}P^2$$

$$P = \sqrt{2\text{m.e.U}}$$

$$\lambda = \frac{h}{\sqrt{2\text{m.e.U}}}$$


#### 1. Electron Diffraction,



$$\tan (2\vartheta) = \frac{D}{2L}$$

$$\frac{\sin 2\theta}{\cos 2\theta} = \frac{D}{2L}$$

$$2 \sin \theta \sim \frac{D}{2L} \text{ for } \theta \to 0$$



$$n\lambda = 2d\sin\theta$$

$$n\lambda = d \times 2 \sin \theta$$

$$\lambda = d \times \frac{D}{2L}; \ n = 1$$

And We have,


$$\lambda = \frac{h}{\sqrt{2\text{m.e.U}}}$$


Now we have relation between D vs U:

$$D = \frac{2.L \cdot h}{d\sqrt{2m \cdot e \cdot U}}$$
L = 13.5 cm

#### 1. Electron Diffraction,

- a) Determination of wavelength of the electrons
- b) Verification of the de Broglie's equation
- c) Determination of lattice plane spacings of graphite





2. Photo Electric effect,

#### 2. Photo Electric effect,

Most of the metals under influence of radiation, emit electrons

Albert Einstein in 1905

- That the emission process depends strongly on frequency of radiation.
- For each metal there exists a critical frequency such that light of lower frequency is unable to liberate electrons, while light of higher frequency always does.
- The emission of electron occurs within a very short time interval after arrival of the radiation and member of electrons is strictly proportional to the intensity of this radiation.
- These facts are strong evident that the energy of the radiation is quantised:

$$E = h\nu$$

h : Plank's constant,  $\nu$  : Frequency

#### 2. Photo Electric effect,

• Energy of the bound electrons in metal is:

$$E = e \cdot \phi$$

 $\phi$ : work function


Then for the emission of the electrons,

$$h\nu > e.\phi$$


With the additional velocity of the electrons,

$$h\nu = \frac{1}{2}mv^2 + e \cdot \phi$$


#### 2. Photo Electric effect,




#### 2. Photo Electric effect,



1-Light source, 2-Guide, 3-Scale, 4-Drawtube, 5-Cover, 6-Focus lens, 7-Vacuum Phototube, 8-Base for holding the Phototube, 9-Digital Meter, 10-Display mode switch, 11-Current multiplier, 12-Light intensity switch, 13-Filter set, 14-Accelerate voltage adjustor, 15-Lens cover, 16-Voltage direction switch, 17-Power switch, 18-Power indicator.



#### 2. Photo Electric effect,



- The light source is used to shine light on a photodiode to generate a photo current.
- We stop the photocurrent by applying a potential with retarding potential technique
- The potential required to stop the photocurrent is called as stopping potential (V<sub>s</sub>), and the kinetic energy of the electrons is defined as  $E_e=\frac{1}{2}mv^2=eV_s$

#### 2. Photo Electric effect,

$$E_e = \frac{1}{2}mv^2 = eV_s$$

Then,

$$h\nu = eV_s + e.\phi$$


$$V_{s} = \frac{h}{e}\nu - \phi$$

$$V_s = -\phi$$
 for  $\nu = 0$ 

Determination of Planck's Constant

Form the slope of the equation:

$$V_{s} = \frac{h}{e}\nu - \phi$$



3. Velocity of Light,

#### 3. Velocity of Light,

A periodic light is an electromagnetic signal which intensity is dependent on time, and change its phase by distance.

$$I = I_0 + \Delta I_0 \cdot \cos(2\pi \cdot \nu \cdot t)$$

We can write simply as alternative signal


$$U = a \cdot \cos(2\pi \cdot \nu \cdot t)$$




#### 3. Velocity of Light,

It will show a phase difference at the receiver as


$$U = a \cdot \cos(2\pi \cdot \nu \cdot t - \Delta\phi)$$






Lets make a path difference for a known phase difference by measuring the phase at the receiver

#### 3. Velocity of Light,



#### 3. Velocity of Light,



The extended path light path is  $\Delta l = 2.\Delta x$ 

For the phase difference ( $\Delta \phi = \pi$ ) the time required is  $\Delta t = 1/2f$ , f is modulation fequency

The velocity of light is 
$$C = \frac{\Delta l}{\Delta t} = 4f. \, \Delta x$$

#### 3. Velocity of Light, Determination of Refractive index:





Fig. 1 In a medium with the refractive index n light propagates at a lower velocity than in vacuum. This leads to a change in the propagation time t of light along a path of length d.

Velocity of light in medium is:

$$C_{M} = \frac{C}{n}$$

3. Velocity of Light, Determination of Refractive index:



#### 3. Velocity of Light, Determination of Refractive index:

Refractive index of the medium is:

$$n = \frac{C}{C_{\rm M}} = \frac{2.\Delta x}{l_m} + 1 + \frac{k.C}{f.l_m}$$

$$\frac{k \cdot C}{f \cdot l_m} \sim 6.k$$
, for 1m water medium

$$\frac{k \cdot C}{f \cdot l_m} \sim 20.k$$
, for 30cm resin medium

For 4. Frank Hertz and 5. Stefan Boltzmann, Look to Prof. Rumi De's notes

#### Todo list in each experiments:

- 1. Electron Diffraction (Estimate the Planck's Constant)
  - 5 Different voltage between 3kV to 5kV
  - For each voltage at least 6 readings of diameter of each ring
- 2. Velocity of Light (Find the speed of light and refractive index of medium)
  - At least 6 Readings of c in air
  - 6 readings in medium
- 3. Photo Electric effect (Estimate the Planck's Constant)
  - For each value of frequency, at least 5 readings averages.
  - Two graphs of stopping potential vs frequency at two different intensities, either by varying distance or bulb current.
- 4. Frank Hertz (Determination of Ar gas first excited state)
  - At least two (three preferable) sets of readings by varying Extraction voltage and retardation voltage
- 5. Stefan Boltzmann (Verify the Stefan Boltzmann Law)
  - Room temperature Resistance of bulb at 6 to 8 different current
  - V vs I
  - At least 6 to 8 reading in bulb glowing condition.

Marks divison (Total 100 Marks)

- 1. Lab notes (25 Marks),
  - Each experiment carried 5M
- 2. Midsem VIVA (25 Marks)
  - After minimum three experiments done by all sub-groups
- 3. Endsem Practical and VIVA (50 Marks),

#### All the best